TorqControl®

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
The success of the implant treatment mostly depends on the precise tightening of the parts, whether surgical or prosthetic, placed directly on implant. A pre-stressed tightening of the screw, will help avoiding any risk of screw loosening while a calibrated tightening will only be guaranteed by the use of a precision instrument, offering a torque control system.

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Functioning principle

The TorqControl is a manual dynamometrical declutching wrench fitted with an adjusting knob allowing a control of the tightening torque from 10 to 35 N.cm.

As soon as the selected torque is reached, that is when the screw is tightened under its recommended tightening value, the wrench frees the mandrel and stops transmitting the tightening effort above the optimised value, tightening being automatically stopped.

Reliability

This wrench has been designed so to guarantee a very high reliability in its use. Indeed, the 7 proposed torque values (10 - 15 - 20 - 25 - 30 - 32 and 35 N.cm), factory calibrated, make it a tool of extreme tightening precision. These values allow moreover the technician and the clinician to proceed with different calibrated tightenings, whether prosthetics in trial phase or pre-stressed, or in association with a good respect of the recommended surgical protocol (cover screw, healing screw).

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PRE-STRESSED ASSEMBLING TECHNIQUE

This technique consists in optimising the tightening effort of the prosthetic part's fixing screw. The tightening torque which will have to be reached to attain this tension is mainly determined according to the elastic limit of the screw.

As such, for a M2 type screw, in Ta6V titanium alloy, with an elasticity limit of circa 1000 MPa (Mega Pascal), the required tightening torque will be evaluated between 32 and 35 N.cm. For a M2 type screw in gold alloy (e.g. Pivozyl® 2), with an elasticity limit of circa 800 MPa, its required tightening torque will be evaluated between 30 and 32 N.cm.

Ergonomy

The torque selection is made very fast directly from the dynamometrical wrench by a simple orientation of the cursor on the adjusting knob in par with the required tightening value. This value remains very easy to modify from one tightening to the next, with a single manipulation.

The access in mouth is optimized thanks to its new micro-head. Its lightness gets a greater facility of use and a better freedom of movement.

Thanks to its length longer than the most traditionnal dynamometrical wrenches, and also due to the optimized orientation at 100° of the mandrels' axis compared with the tool's axis, the ergonomic of the Torq Control® has also been thought so to favor the access to the screws placed in the posterior areas.

Finally, its design made in one piece with large smooth surfaces eases its maintenance, during the cleaning phases as well as when sterilised in autoclave at 135°C.

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


■ Applications

Tightening of the cover screw

- Diameter M2
- Recommended torque : 10 N.cm*
- Allows a secure tightening in the difficult or less difficult areas, as well as an easier withdrawal of the screw when entering a second surgical time while reducing any risk of cold welding of the screw on the implant.

Tightening of the healing screw

- Diameter M2
- Recommended torque : 10 N.cm*
- Allows the control of the torque indirectly transmitted by the screw on the implant during its placing (especially during 1 surgical time surgeries) in order to preserve the primary stability of the implant. Also allows a secure tightening in the difficult or less difficult areas, as well as an easier withdrawal of the screw while reducing any risk of cold welding of the screw on the implant.

Tightening of direct abutment on implant, titanium abutment, overcast abutment.

- Diameter M2 gold or titanium
 - Recommended torque of 30 to 35 N.cm*
- These values are depending upon the size of the thread's pitch, the number of threads, the type of gold and titanium alloy. Please contact the manufacturer of the screws.

Tightening of the gingival abutment

- Ø M2 gold or titanium
 - recommended torque of 30 to 35 N.cm*
- These values are depending upon the size of the thread's pitch, the number of threads, the type of gold and titanium alloy. Please contact the manufacturer of the screws.

Tightening of prosthesis on abutment

- Diameter M1.4 Gold or Titanium
- Recommended torque : 15 N.cm*
- Precision and safety of the tightening.

TorqControl®

anthogyr