NiTiControl®

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
*The use of the Nickel-Titanium files is nowadays common practice in terms of endodontia. Thanks to an easy-to-use protocol and to an extreme flexibility of the file, the canal route is perfectly respected. Suited to these new files, the technical characteristics of the NiTi Control allow a **better monitoring of the continuous rotation** and a control of the possible parietal restraints. Its **automatic declutching** function represents **additional safety**.*

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

Thanks to a new chrome plated coating and to an even surface, it is easier to clean the NiTi Control® contra-angle. It is also proved to be more resistant to repeated decontaminations and sterilisations. The life span of the contra-angle is extended without altering neither its lightness nor its handiness.

Working Method

• Continuous rotation

The NiTi Control contra-angle is aimed at endodontic preparations using the continuous rotation technique. The NiTi Control contra-angle does not require any specific motor: it is directly connected to the Unit.

It is available in two versions :

- 64:1 for the 20 000 rpm motor (maximum speed)
- 128:1 for the 40 000 rpm motor (maximum speed).

• Adjustable torque and automatic declutching

The NiTi Control contra-angle is equipped with a torque adjusting ring located at head level. Each adjustment relates to a different torque value. According to the diameter and the length of the file, the most adapted adjustment is used, in other words, the maximum torque which could be sustained by the file is used. As soon as the value is reached, the rotation is automatically stopped.

Advantages

Maximum visibility

→ Its micro push button head enables an easy access and its dimensions are ideal for the use of the «mini» files

Safety

→ The canal treatment is made safer due to the declutching function

Faster

→ Insertion of the file without operating the push button
→ Ease of the push button.