

EOFY PROMOTION

MAY - JUNE 2024



PRIME LINE



AeroAcoustic blade



Glass rod fiber optics

25 W

Powerful & Precise Power Head - 25W



Reverse flow protection



Automatic pressure regulation¹

55 dB

Incredibly Quiet Only 55 dB noise



Refined ceramic bearings



Scratch resistance
Iconic Prime body



Available for
KaVo*, NSK*, W&H*
and Sirona* Couplings



**BUY 2
&
GET 1
FREE**

\$1,683.00



**Now with 2 Years
FairCare Warranty²**

Unconditional warranty
except physical damage

Coupling Connection

Kavo
NSK
W&H
Sirona

Power Head 25 Watt

HP21KL
HP21NL
HP21WL
HP21SL

Small Head 18 Watt

HP22KL
HP22NL
HP22WL
- N/A -

COUPLINGS



✓ MK-dent & Kavo with LED & water
• **QC6016KW**
\$ 495.00



✓ NSK connection with LED
• **QC6016NT**
\$ 385.00



✓ W&H connection with fibre optic
• **QC5016W**
\$ 385.00



✓ Sirona connection with LED & water
• **QC6016SW**
\$ 495.00

✓ Non-optic
• **QC4014K**
\$ 296.00

✓ Non-optic
• **QC4014NT**
\$ 296.00

✓ Non-optic
• **QC4014W**
\$ 296.00

✓ Non-optic
-NA-

1) Only available for Kavo coupling 2) FairCare warranty applicable on Prime Line high speeds sold after 1st Oct 2018

**Buy 3 - Get 20% OFF
Mix & Match**


**NEW
RANGE**

**With 2 Years
FairCare Warranty**





**SAFE-TECH
PUSH
BUTTON**



 1:5 Red ●


 Max. 200,000 rpm

 Internal water

 For FG-burs
(Ø 1.6 mm)


● LP15L (F/Optic)
RRP: \$ 1,826.00

Buy 3 - \$1,460.80 ea

 1:1 Blue ●

 Max. 40,000 rpm


 Internal water


 For CA-burs
(Ø 2.35 mm)


● LP11L (F/Optic)
RRP: \$ 1,287.00

Buy 3 - \$1,029.60 ea

 5.4:1 Green ●

 Max. 7,407 rpm

 Internal water


 For CA-burs
(Ø 2.35 mm)


● LP54L (F/Optic)
RRP: \$ 1,320.00

Buy 3 - \$1,056.00 ea

 1:1 Blue ●

 Max. 40,000 rpm

 Internal water

 For HP-burs
(Ø 2.35 mm)

● LP01L (F/Optic)
RRP: \$ 1,045.00

Buy 3 - \$836.00 ea

AIR MOTORS 



NEW
Shorter
and
Lighter

- ✓ Internal water
- ✓ LED

● AM1006
\$ 990.00



- ✓ Internal water

● AM1004
\$ 748.00



- ✓ External water

● AM0014
\$ 473.00



● AM0004
\$ 473.00

NEW

KAVO COMPATIBLE HEADS 



- Kavo 68LH compatible
- Blue (1:1)
- LCH11
\$ 473.00



- Kavo 67LH compatible
- Green (2:1)
- LCH21
\$ 473.00



- Kavo 31LR compatible
- Green (2:1)
- LCH21P
\$ 484.00

ECO LINE TURBINES



Now with 1 Year FairCare Warranty¹
Unconditional warranty except physical damage

Now with
24W²



- 62 dB Operating Noise
- AeroPerformance blade
- Refined ceramic bearings
- Glass rod fiber optics
- Smooth Eco body



Available for
KaVo*, NSK*, W&H*
and Sirona* Couplings

Savings from \$561



BUY 3
or MORE
Only
\$792
each

3 x Eco Line HS
with F/Optic
for any
connection

Coupling Connection	Power Head 24W (F/Optic)	Small Head 16W (F/Optic)	Power Head 24W (N/Optic)	Small Head 16W (N/Optic)
KaVo	HE21KL	HE22KL	HE21K	HE22K
NSK	HE21NL	HE22NL	HE21N	HE22N
W&H	HE21WL	HE22WL	HE21W	HE22W
Sirona	HE21SL	- N/A -	- N/A -	- N/A -

• \$979.00 (F/Optic)

• \$693.00 (N/Optic)

Buy 3 - Get 10% OFF

ECO LINE POWER EDITION



NEW

**28W
POWER EDITION**

- 28 Watt Power
- 60 dB Operating Noise
- AeroPerformance blade
- Glass rod fiber optics
- Refined ceramic bearings
- Smooth Eco body

• HE20KL
\$1,265.00

Buy 3 or more
Only \$1,089 ea



With 18 months
FairCare Warranty
Unconditional warranty
except physical damage

Available for only
KaVo Multiflex coupling fitting

ECO LINE FIXED CONNECTION TURBINES

- Direct 4 hole connection
- Refine ceramic bearings
- Push Button System
- HE17 - Standard head with 4 hole water spray
- HE18 - Small head with 3 hole water spray



\$539.00

Buy 3 - Get 10% OFF

1) FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019 2) Eco Line Power Head 24 Watt

Buy 3 - Get 1 FREE

- Mix & Match
- Cheapest Free



FairCare Warranty¹
 1 year unconditional warranty
 except physical damage



⚙️ 1:5 Red ●
 ○ Internal water
 For FG-burs
 (∅ 1.6 mm)
 ● LE15L (F/Optic)
 RRP: \$1,595.00

Buy 3 - Get 1 FREE

● LE15 (Non-optic)
 RRP: \$1,375.00
Buy 3 - Get 1 FREE

⚙️ 1:1 Blue ●
 ○ Internal water
 For CA-burs
 (∅ 2.35 mm)
 ● LE11L (F/Optic)
 RRP: \$1,155.00

Buy 3 - Get 1 FREE

● LE11 (Non-optic)
 RRP: \$891.00
Buy 3 - Get 1 FREE

⚙️ 5.4:1 Green ●
 ○ Internal water
 Screw-in & Snap-on
 Prophy Cups
 ● LE54LP (F/Optic)
 RRP: \$1,287.00

Buy 3 - Get 1 FREE

● LE54P (Non-optic)
 RRP: \$1,045.00
Buy 3 - Get 1 FREE

⚙️ 1:1 Blue ●
 ○ Internal water
 For HP-burs
 (∅ 2.35 mm)
 ● LE01L (F/Optic)
 RRP: \$990.00

Buy 3 - Get 1 FREE

● LE01 (Non-optic)
 RRP: \$825.00
Buy 3 - Get 1 FREE

BASIC LINE CONTRA - ANGLES 



⚙️ 1:1 Blue ●
 ○ Push button
 ● LB11/11
 Buy 1 - \$418.00

Buy 3: \$339.00 ea
Buy 4+ \$319.00 ea

⚙️ 1:1 Blue ●
 ○ Latch Type
 ● LB11/11C
 Buy 1 - \$418.00

Buy 3: \$339.00 ea
Buy 4+ \$319.00 ea

⚙️ 4:1 Green ●
 ○ Snap-on Cups
 ● LB41/11P
 Buy 1 - \$594.00

Buy 3: \$484.00 ea
Buy 4+ \$473.00 ea

⚙️ 1:1 Blue ●
 ○ For HP-burs
 ● LB02
 Buy 1 - \$418.00

Buy 3: \$339.00 ea
Buy 4+ \$319.00 ea

1) FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019

✓ ULTRASONIC SCALER

\$ Buy 3 - Get 1 FREE

1 Year Warranty



EMS Non-optic Type EN-041

- SC21E
- Buy 1 - \$462.00
- Buy 3 - Get 1 FREE**



EMS Piezon Non-optic Type EN-061

- SC21E2
- Buy 1 - \$561.00
- Buy 3 - Get 1 FREE**



EMS Piezon LED Type EN-060A

- SC21EL
- Buy 1 - \$748.00
- Buy 3 - Get 1 FREE**



KaVo Piezo LED Type 1.007.3995

- SC21KL* **NEW**
- Buy 1 - \$748.00
- Buy 3 - Get 1 FREE**



Acteon-Satelec Non-optic Type F 12 281

- SC21S
- Buy 1 - \$462.00
- Buy 3 - Get 1 FREE**



Acteon-Satelec LED Type F 12 609

- SC21SL
- Buy 1 - \$748.00
- Buy 3 - Get 1 FREE**

*Only for use with EMS compatible scaler tips

✓ SCALER TIPS



Buy 3 - Get 1 FREE while stock lasts

- ✓ High-quality Diamond-Like-Carbon coating
- ✓ Extremely hard and wear-resistant surface

• **SCALING DLC COATED** Only available for Satelec fitting

ST03SD
General Scaling
\$88.00

ST04SD
SATALEC 2
\$132.00

• **PERIO DLC COATED** Only available for Satelec fitting

ST06SD
SATELEC 10Z
\$132.00

ST10SD
SATALEC 10P
\$154.00

• **ENDO DLC**

ST08SD
SATALEC S07
\$154.00

• **SCALING STAINLESS STEEL** Only available for Satelec fitting

ST02S
SATELEC 1S
\$63.80

ST03S
General Scaling
\$63.80

ST04S
SATELEC 2
\$88.00

• **PERIO STAINLESS STEEL** Available for EMS & Satelec fittings

ST07E/ST07S
EMS DPL3
SATELEC 10X
\$110.00

ST09E/ST09S
EMS HPL
\$110.00

ST10E/ST10S
SATALEC 10P
\$110.00

ST11E/ST11S
EMS RT1
\$110.00

Universal Torque Wrench



TORQ CONTROL®

Precision instrument with a Torque Adjustment System for tightening of Prosthetic components at specific Torque, for any implant system.



Automatic Declutching when Selected Torque is reached

- Compliance with Protocols for screw tightening
- Limits the risk of the screw fracturing or loosening




Universal connection of Prosthetic mandrels

- 7 different Torque levels (from 10 to 35 N.cm)
- Can be used with ALL Implant Systems



Micro-Head and Ergonomic Handle

- Easy access in restricted area with Excellent Visibility



Easy and Quick to use

- Reduced motion compared with using a ratchet wrench while tightening



\$1,210.00


Lightweight
ONLY 135 gms

Order code: 15501

3 Easy Steps Process



1. Please select the recommended Torque



2. Start rotating



3. Screw automatically stops when selected torque is reached



www.anthogyr.com

Anthogyr
A Straumann Group Brand

▼ PROPHY LINE AIR POLISHING DEVICE



Buy 1 Polisher and receive
• 1 bottle of Propphy Powder FREE
and
• 50% discount on 1 additional tip

• PR1011
\$1,375.00



• Ergonomically shaped, compact and light-weight air-driven polishing device

• 360° rotatable outer handle due to unique RotaGrip® for simple and fatigue-proof operation



• Propphy Powder (300 gms)
\$ 60.50 each

• Lemon - PP1012

Additional Propphy tip



• SP3011
\$ 329.00

Coupling connection order codes:

Kavo - PR1011K, NSK - PR1011N, W&H - PR1011W, Sirona - PR1011S

▼ LEDs

- ✓ Quality illumination, No Glare
- ✓ 5500K Natural Daylight Colour
- ✓ Clear Vision & Sharper Contrast
- ✓ Longer Durability & Economical



\$88.00

- BU8012 - KaVo & MK-dent coupling
- BU8012KM - KaVo & MK-dent motors
- BU8012WH - W&H handpieces
- BU8012WLT - W&H airmotors LT25
- BU8012WEA - W&H motors 40LT, 50 LT & A51LT
- BU8012SG - Sirona coupling and motors
- BU8012N - NSK coupling
- BU8012NT - NSK motors TiM40 & Ti205L
- BU8012NM - NSK coupling KCL,SCL,PTL LED

▼ PREMIUM LUBRICATION OIL



- ✓ Fully Synthetic
- ✓ For any brand of handpieces and motors

• LU1011 (500 ml)
\$ 68.20



- ✓ Fully Synthetic
- ✓ For KaVo Quattrocare 2 lubrication system

• LU1022 (500 ml)
\$ 68.20



- ✓ For W&H Assistina lubrication system

• LU1031 (500 ml)
\$ 88.00

MANUFACTURER



MK-dent GmbH
 Marie-Curie-Str. 2
 22941 Bargteheide
 GERMANY
www.mk-dent.com

YOUR SPECIALIST DEALER

• All catalogue prices are inclusive of GST • Promotion valid until 12th April 2024 or while stocks last

KaVo, W&H, Sirona, NSK, Bien-Air, EMS and Acteon-Satelec are registered trademarks. MK-dent GmbH has no economic connection with the aforementioned companies.