


STAINLESS STEEL SYRINGES

STERILIFE®



- 100% stainless steel.
- Highly resistant to sterilisation.
- Stainless steel syringes with lateral insertion of the cartridge.



Ref. 2310	Syringe without aspiration, 1,8 ml Metric thread
Ref. 2310US	Syringe without aspiration, 1,8 ml American thread
Ref. 2310-22US	Syringe without aspiration, 2,2 ml American thread
Ref. 2320	Syringe with self-aspiration, 1,8 ml Metric thread
Ref. 2320US	Syringe with self-aspiration, 1,8 ml American thread
Ref. 2320-22US	Syringe with self-aspiration, 2,2 ml American thread
Ref. 2330	Syringe with manual aspiration, 1,8 ml Metric thread
Ref. 2330US	Syringe with manual aspiration, 1,8 ml American thread
Ref. 2330-22US	Syringe with manual aspiration, 2,2 ml American thread




CHROME PLATED SYRINGES

 → The syringes with aspiration are fitted with a cork screw.


Ref.1640	Straight handle syringe without aspiration , 1,8 ml Metric thread
Ref.1640US	Straight handle syringe without aspiration, 1,8 ml American thread
Ref. 3800	Straight handle syringe with aspiration, 1,8 ml Metric thread
Ref. 3800US	Straight handle syringe with aspiration, 1,8 ml American thread
Ref. 3800A	Thumb ring handle syringe with aspiration, 1,8 ml Metric thread
Ref. 3800AUS	Thumb ring handle syringe with aspiration, 1,8 ml American thread

Accessories

Ref.3804F	Corkscrew for aspiration
-----------	--------------------------



INTRALIGAMENTAL SYRINGES





Ergonomically designed, these technical syringes allow every kind of anaesthesia, especially intraligamental and intraseptal.

- Progressive mechanism which makes the injection very accurate and smooth.
- A system of levers is used to reduce the force required to administer the injection:
 - Ergoject x3.
 - Miniject x5.
- Two tips available : straight tip and 30° angulated tip with needle insert cone.
- Supplied with 2 disposable protective tubes.

Ref. 2300	Miniject syringe with angulated barrel, 1,8 ml metric thread
Ref. 2300US	Miniject syringe with angulated barrel, 1,8 ml American thread
Ref. 2400	Ergoject syringe with straight barrel, 1,8 ml metric thread
Ref. 2400US	Ergoject syringe with straight barrel, 1,8 ml American thread

Accessories

Ref. 2390	Straight barrel	
Ref. 2392	Angulated barrel	
Ref. 2429	Pack of 12 protective plastic tubes	