

PRIME LINE - The Premium Range

Buy 3 for \$946.⁶⁷ each **SAVE \$1,350.⁰⁰**
Buy 2 for \$1,047.⁰⁰ each **SAVE \$700.⁰⁰**



HP22KL



HP21KL



Power-Torque head **25W**
and Small head **18W**



Ceramic bearings for longer
life and quiet operation



Premium spindle with carbide
insert for high concentricity



40,000 LUX solid fibre
optic cellular glass rod



Pressure reduction system
in KaVo^{*} Multiflex fitting



Scratch resistance
titanium coated prime-grip



Integrated non-return valve
to prevent cross-infection



Efficient cooling with
4-port or 3-port spray

\$1,397.00 each

Coupling connection	Power Head	Small Head
KaVo [*]	HP21KL	HP22KL
NSK [*]	HP21NL	HP22NL
W&H [*]	HP21WL	HP22WL
Sirona [*]	HP21SL	- NA -

German products
at affordable prices
and backed by
factory warranty

All prices in this catalogue are inclusive of GST

All Specials/Offers valid until
30th September 2018

www.mk-dent.com

^{*} KaVo, Sirona, NSK, and W&H are registered trademarks of KaVo Dental GmbH, Sirona Dental Systems GmbH, NSK, and W&H Dental Work Bürmoos Inc. respectively. MK-dent is not affiliated in any way with these companies.



FREE

DLC coated tip
(similar to EMS-A type)
worth \$110
with
every scaler

EMS® connection

~ **SC21E (SETSC2)** **\$385.00**



EMS® N/Optic
connection type EN-041

~ **SC21E2 (SETSC5)** **\$495.00**



EMS® Piezon® NO-PAIN
N/Optic connection EN-061

~ **SC21EL (SETSC1)** **\$594.00**



EMS® Piezon® LED NO-PAIN
connection type EN-060A

23 different types of tips available
for EMS® & Satelec® Scalers
Please visit www.mk-dent.com

Satelec® connection

~ **SC21S (SETSC4)** **\$385.00**



Acteon-Satelec® N/Optic
connection type F 12 281

~ **SC21SL (SETSC3)** **\$594.00**



Acteon-Satelec® LED
connection type F 12 609

KaVo® connection

~ **SC21KL* (SETSC6)** **\$594.00**



KaVo® PiezoLED®
connection Ref 1.007.3995

*Only for use with MK-dent STxxE/ED
or EMS® original scaler tips



SAFE-TECH
Push button

NEW

- **LE15L (Optic)**
\$1,292.50
- Speed increasing (1:5)
- Max 200,000 rpm
- Ceramic bearings
- Safe-Tech Push button
- Internal water
- For FG-burs (Ø 1.6mm)

SAVE \$260

Buy 2 for \$1,162.50 ea

ECO LINE

SAVE \$630

Buy a set of 3
for \$3,245.00

That's only
\$1,081.65 each



40,000 LUX solid fibre
optic cellular glass rod



Power head **21W** and
Small head **16W** options



Ceramic bearings for longer
life and quiet operation



Efficient cooling with
4-port or 3-port spray



Smooth grip with modern
and ergonomic design



ECO LINE

SAVE \$250

Buy a set of
3 for
\$1,895⁰⁰

**That's only
\$631.65 each**

Buy one for \$715.00

Coupling connection	Power Head	Small Head
KaVo*	HE21KL	HE22KL
NSK*	HE21NL	HE22NL
W&H*	HE21WL	HE22WL
Sirona*	HE21SL	- NA -

ECO LINE

**Buy 2 and Get 10% OFF
Mix & Match**



- **LE01L (Optic)**
\$825.00
- **LE01 (N/Optic)**
\$660.00

- Blue (1:1)
- Max 40,000 rpm
- Internal water
- For HP-burs
(Ø 2.35mm)



- **LE11L (Optic)**
\$825.00
- **LE11 (N/Optic)**
\$660.00

- Blue (1:1)
- Max 40,000 rpm
- Internal water
- For CA-burs
(Ø 2.35mm)



- **LE54PL (Optic)**
\$825.00
- **LE54P (N/Optic)**
\$660.00

- Green (5.4:1)
- Max 7,407 rpm
- Internal Water
- For Snap-on &
Screw-in Prophy
attachments

NEW

UPGRADED
MODEL

**Buy 4
Get 1 FREE**
Mix & Match

BASIC LINE



LB11/11
\$385.00

- 20,000 rpm
- Blue (1:1)
- Push-button Head



LB41/11P
\$385.00

- 5,000 rpm
- Green (4:1)
- For Snap-On Prophylaxis Cups



LB41/11
\$440.00

- 5,000 rpm
- Green (4:1)
- Push-button Head



LB01
\$385.00

- 20,000 rpm
- Blue (1:1)

AIR MOTORS



AM1116
\$990.00

- Internal water
- With LED



AM1004
\$605.00

- Internal water
- Non-optic



AM0004
\$385.00

- No water
- Non-optic



AM0014
\$385.00

- Ext. water
- Non-optic

COUPLINGS

QC6016K
\$440.00



for KaVo* &
MK-dent
with LED

QC6016KW
\$495.00



for KaVo* &
MK-dent
LED & water

QC6016NT
\$330.00



for NSK*
connection
with LED

QC6016SW
\$440.00



for SIRONA*
connection
LED & water

QC5016W
\$330.00



for W&H*
connection
with optic