

PRIME LINE
Premium Handpieces



Incredibly Quiet !!
55dB Operating Noise



Powerful & Precise !!
Power-Torque Head 25W



NEW LONG LIFE
Premium spindle with carbide insert for high concentricity



Pressure reduction system in KaVo® Multiplex fitting



Integrated non-return valve to prevent cross-infection



Ceramic bearings for longer life and quiet operation



40,000 LUX solid fibre optic cellular glass rod



Scratch resistance
Titanium coated Prime-grip



Efficient cooling with
4-port or 3-port spray



UPGRADED MODEL

HP21NL

HP21KL



FairCare
Cost Effective and Comprehensive Care for Peace of Mind

- Comprehensive & Unconditional 2 Year FairCare* warranty including parts and labour
- Affordable cost of genuine spare parts for post warranty

Available coupling connection	Power Head 25W	Small Head 18W
KaVo®	HP21KL	HP22KL
NSK®	HP21NL	HP22NL
W&H®	HP21WL	HP22WL
Sirona®	HP21SL	- NA -

\$1,397.00 each

- Buy 3 - \$995.⁵⁰ ea | **SAVE \$1,204.⁵⁰**
- Buy 2 - \$1,067.⁰⁰ ea | **SAVE \$660.⁰⁰**

* 2 years FairCare unconditional warranty (except physical damage) applicable on Prime Line high speeds sold after 1st Oct 2018

KaVo, Sirona, NSK, and W&H are registered trademarks of KaVo Dental GmbH, Sirona Dental Systems GmbH, NSK, and W&H Dental Work Bürmoos Inc. respectively. MK-dent is not affiliated in any way with these companies.



Gulmohar Dental Pty Ltd
14/1 Ricketts Road, Mt Waverley, VIC 3149
Ph: 03 9544 6980 Fax: 03 9544 6946
www.gulmohar.com.au

- All catalogue prices are inclusive of GST
- Specials / Offers valid until 30th April 2019

*In relation to any promotions or offers by, relating to or referencing Credabl Pty Ltd (ACN 615 968 100) (Credabl) that are contained in this document (Credabl Promotion), this document is a guide only and does not constitute any recommendation on behalf of Credabl. The information in this document is general in nature and Credabl has not taken into account your personal objectives or financial circumstances or needs when preparing it. Before acting on this information you should consider if it is suitable for your personal circumstances. Credabl is not offering financial, tax or legal advice. You should obtain independent financial, tax or legal advice as appropriate. Credabl reserves the right to cease offering these products at any time without notice. The Credabl Promotion is subject to Credabl's terms and conditions, fees and charges and eligibility criteria. MK-Dent Australia Pty Ltd **ACN 161 945 946 (MK-Dent)** is not an agent or partner of Credabl and is a separate and independent entity from Credabl. In no event will Credabl be liable for any actions, representations and / or decisions of MK-dent or its related bodies corporate representatives, directors, officers, employees, agents or advisors. **The issuer and credit provider of the Credabl Promotion is Credabl Pty Ltd (ACN 615 968 100) Australian Credit Licence No. 499547.**

ULTRASONIC SCALER

MK-dent
GERMANY · SINCE 1995



**ONE YEAR
WARRANTY**

**Buy 3 - Get 15% OFF
PLUS
FREE tip (worth \$110.00)
with every scaler handpiece**

EMS® connection



EMS® Piezon® LED No Pain
connection type EN-060A

SC21EL
(SetSC1)

- Buy 1 \$594.00
- **Buy 3 \$504.90 ea**



EMS® Piezon® No Pain Non-Optic
connection type EN-061

SC21E2
(SetSC5)

- Buy 1 \$495.00
- **Buy 3 \$420.75 ea**



EMS® Non-Optic
connection type EN-041

SC21E
(SetSC2)

- Buy 1 \$385.00
- **Buy 3 \$327.25 ea**

KaVo® connection

Coming soon !



KaVo® PiezoLED® connection
*Only for use with MK-dent STxxE/ED or
EMS® original scaler tips

SC21KL*
(SetSC6)

- Buy 1 \$594.00
- **Buy 3 \$504.90 ea**

Acteon®/Satelec® connection



Acteon-Satelec® Non-Optic
connection type F 12 281

SC21S
(SetSC4)

- Buy 1 \$385.00
- **Buy 3 \$327.25 ea**



Acteon-Satelec® with LED
connection type F 12 609

SC21SL
(SetSC3)

- Buy 1 \$594.00
- **Buy 3 \$504.90 ea**

































SC21KEY
Torque Key
\$38.50

NEW

**Innovative DLC coated Black tips
for longer life and optimum visibility**

**Buy 4 - Get 10% OFF
Mix & Match**

	 Satelec® 1 \$77.00	 Satelec® 1 \$110.00	 EMS® P Satelec® 1S \$77.00	 EMS® P Satelec® 1S \$110.00	 EMS® A \$77.00	 EMS® A \$110.00
EMS® / Satelec®	ST01E / ST01S	ST01ED / ST01SD	ST02E / ST02S	ST02ED / ST02SD	ST03E / ST03S	ST03ED / ST03SD
Application	Scaling	Scaling	Scaling	Scaling	Scaling	Scaling
Coating	Stainless Steel	DLC	Stainless Steel	DLC	Stainless Steel	DLC
	 EMS® C Satelec® 2 \$77.00	 EMS® C Satelec® 2 \$110.00	 EMS® D Satelec® 5AE \$77.00	 EMS® PS Satelec® 10Z \$88.00	 EMS® PS Satelec® 10Z \$121.00	 EMS® DPL3 Satelec® 10X \$99.00
EMS® / Satelec®	ST04E / ST04S	ST04ED / ST04SD	ST05E / ST05S	ST06E / ST06S	ST06ED / ST06SD	ST07E / ST07S
Application	Scaling	Scaling	Endodontics	Periodontal	Periodontal	Periodontal
Coating	Stainless Steel	DLC	Stainless Steel	Stainless Steel	DLC	Stainless Steel
	 Satelec® S07 \$99.00	 Satelec® S07 \$143.00	 EMS® HPL \$99.00	 Satelec® 10P \$99.00	 Satelec® 10P \$132.00	 EMS® RT1 \$99.00
EMS® / Satelec®	ST08E / ST08S	ST08ED / ST08SD	ST09E / ST09S	ST10E / ST10S	ST10ED / ST10SD	ST11E / ST11S
Application	Periodontal	Periodontal	Periodontal	Periodontal	Periodontal	Periodontal
Coating	Stainless Steel	DLC	Stainless Steel	Stainless Steel	DLC	Stainless Steel
	 120° File Holder \$88.00	 90° File Holder \$88.00	 180° File Holder \$88.00	 Satelec® ET40 \$99.00	 Satelec® ET40D \$99.00	 Satelec® ET20 \$99.00
EMS® / Satelec®	ST12E / ST12S	ST13E / ST13S	ST14E / ST14S	ST15E / ST15S	ST16E / ST16S	ST17E / ST17S
Application	Endodontics	Endodontics	Endodontics	Endodontics	Endodontics	Periodontal
Coating	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
	 Satelec® ET20D \$99.00	 EMS® SB Satelec® EX1 \$110.00	 EMS® SBM Satelec® EX2 \$110.00	 EMS® SBD Satelec® EX3 \$110.00	 Satelec® EXL \$110.00	 Satelec® EXR \$110.00
EMS® / Satelec®	ST18E / ST18S	ST19E / ST19S	ST20E / ST20S	ST21E / ST21S	ST22E / ST22S	ST23E / ST23S
Application	Cavity Preparation	Cavity Preparation	Cavity Preparation	Cavity Preparation	Cavity Preparation	Cavity Preparation
Coating	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel

EMS®, Acteon® and Satelec® are registered trademarks of their respective companies. MK-dent is in no way affiliated with EMS® and Acteon®/Satelec®



40,000 LUX solid fibre
optic cellular glass rod



Power-Torque head **21W**
and Small head **16W**



Ceramic bearings for longer
life and quiet operation



Efficient cooling with
4-port or 3-port spray



Smooth-grip with modern
ergonomic design



Made in
Germany



LED
Compatible
Handpieces

ECO LINE

\$715.00 each

Buy 3 - \$632.⁵⁰ ea

Available
coupling
connection

with F/Optic

KaVo®
NSK®
W&H®
Sirona®

Power
Head

Small
Head

HE21KL
HE21NL
HE21WL
HE21SL

HE22KL
HE22NL
HE22WL
- NA -



Now with
1 YEAR

FairCare*

**UNCONDITIONAL
WARRANTY**

NEW



SAFE-TECH
Push button

• LE15L (Optic)

Buy 1 \$1,292.50

Buy 3 \$1,098.60 ea

• LE15 (Non-Optic)

Buy 1 \$1,155.00

Buy 3 \$ 981.75 ea

- Speed increasing (1:5)
- Max 200,000 rpm
- Ceramic bearings
- Safe-Tech Push button
- Internal water
- For FG-burs (Ø 1.6mm)

Buy 3 - Get 15% OFF
Mix & Match



• LE11L (Optic)

Buy 1 \$825.00

Buy 3 \$701.25 ea

• LE11 (Non-Optic)

Buy 1 \$660.00

Buy 3 \$561.00 ea

- Blue (1:1)
- Max 40,000 rpm
- Internal water
- For CA-burs (Ø 2.35mm)



• LE01L (Optic)

Buy 1 \$825.00

Buy 3 \$701.25 ea

• LE01 (Non-Optic)

Buy 1 \$660.00

Buy 3 \$561.00 ea

- Blue (1:1)
- Max 40,000 rpm
- Internal water
- For HP-burs (Ø 2.35mm)

* 1 year FairCare unconditional warranty (except physical damage) applicable on Eco Line high & low speeds sold after 1st Jan 2019

PROPHY LINE

- Ergonomically shaped, compact and light-weight Air-driven polishing device
- 360° rotatable outer handle due to unique RotaGrip® for simple & fatigue-proof operation



PR1011

\$1,155.00

Available coupling connections:

KaVo® - PR1011K W&H® - PR1011W
Sirona® - PR1011S NSK® - PR1011N
BienAir® - PR1011B

• One Year Standard Warranty

FREE
3 bottles of
Prophy Powder
with every
Prophy Line Unit



PP1011
Mint



PP1012
Lemon



PP1013
Strawberry

Prophy Powder

\$60.50 each

PREMIUM OIL

- Fully Synthetic Premium Service Oil •



LU1011
500ml

For any
brand of airmotors
and handpieces

\$55.00



LU1022
500ml

For KaVo®
QUATTROcare 2
lubrication system

\$55.00



LU1031
500ml

For
W&H®
ASSISTINA

\$82.50

BASIC LINE

UPGRADED
MODEL



LB11/11
\$385.00

- 20,000 rpm
- Blue (1:1)
- Push-button Head



LB41/11P
\$385.00

- 5,000 rpm
- Green (4:1)
- For Snap-On Prophy Cups



LB41/11
\$440.00

- 5,000 rpm
- Green (4:1)
- Push-button Head



LB01
\$385.00

- 20,000 rpm
- Blue (1:1)

One Year
Standard
Warranty

Buy 4 - Get 1 FREE
Mix & Match (Cheapest Free)



Made in
Germany

AIR MOTORS

NEW



AM1116
\$990.00

- Internal water
- Titanium Coating
- With LED



AM1004
\$605.00

- Internal water
- Non-optic



AM0004
\$385.00

- No water
- Non-optic



AM0014
\$385.00

- External water
- Non-optic

One Year
Standard
Warranty

COUPLINGS

QC6016K
\$440.00



for KaVo® &
MK-dent
with LED

QC6016KW
\$495.00



for KaVo® &
MK-dent with
LED & water

QC6016NT
\$330.00



for NSK®
connection
with LED

QC6016SW
\$440.00




for SIRONA®
connection with
LED & water

QC5016W
\$330.00



for W&H®
connection
with optic

3 easy steps



Universal Torque Wrench

TORQ CONTROL®


Precision instrument with adjustable torque settings to place prosthetic components at a specific torque

Suitable for any
implant system



Made in France
by

Anthogyr
PRIME MOVER IN ORTHOPAEDY
www.anthogyr.com



Automatic declutching when selected torque is reached

- Compliance with screw tightening protocols
- Limits the risk of screw fracturing or loosening



Universal connection of prosthetic mandrels

- 7 different Torque levels - from 10 to 35 N.cm
- Can be used with ALL Implant systems



Micro-Head and Ergonomic handle

- Excellent visibility and easy access in restricted areas



Easy and Quick to use

- Reduced motion compared with using a ratchet wrench while tightening