

DISCOVER

MK-dent[®]
GERMANY



ECO LINE

HIGH & LOW SPEED HANDPIECES

www.MK-dent.com

ECO LINE



Now also with Multiflex® Connection:
High Flexibility, low Price!

- Smooth Grip
- Chrome coating
- Pushbutton system

- Spindle with steel-insert
- Ceramic bearings
- Glass rod light conductor

Available with Multiflex® connection:

	HE21KL	HE21K	HE22KL	HE22K
Head size	Standard	Standard	Small	Small
Light	✓	—	✓	—
Power (Watt)	18 W	18 W	16 W	16 W
Spray system	4-hole	4-hole	3-hole	3-hole
Connection	Multiflex®	Multiflex®	Multiflex®	Multiflex®

Available with tubing connection:

	HE11	HE12	HE14	HE15L	HE16L	HE17	HE18
Head size	Standard	Standard	Mini	Standard	Small	Standard	Small
Light	—	—	—	✓	✓	—	—
Power (Watt)	18W	18 W	15 W	18W	16W	18W	16W
Spray system	1-hole	1-hole	1-hole	4-hole	1-hole	4-hole	3-hole
Connection	2-hole	4-hole	4-hole	5-hole	5-hole	4-hole	4-hole

ECO LS LINE



Available in the following versions:

	LE11L LE11	LE16L LE16	LE21L LE21	LE27L LE27
Light	with light without light	with light without light	with light without light	with light without light
Transmission	1:1	1:1	2:1	2,7:1
Label	blue	blue	green	green
Spray system	internal	internal	internal	internal
Max. speed	40.000 rpm	40.000 rpm	20.000 rpm	14.814 rpm
Connection	INTRAmatic®	INTRAmatic®	INTRAmatic®	INTRAmatic®
Suitable burs	Contra angle burs	FG burs	Contra angle burs	Contra angle burs

	LE54L LE54	LE74L LE74	LE148L LE148
Light	with light without light	with light without light	with light without light
Transmission	5,4:1	7,4:1	14,8:1
Label	green	green	green
Spray system	internal	internal	internal
Max. speed	7.407 rpm	5.406 rpm	2.703 rpm
Connection	INTRAmatic®	INTRAmatic®	INTRAmatic®
Suitable burs	Contra angle burs	Contra angle burs	Contra angle burs

MK-dent[®]
G E R M A N Y



MK-dent[®]
G E R M A N Y

Eichenweg 7b
D-22941 Bargteheide, Germany

Phone: 0049 (0)4532 / 400 49-0

Fax: 0049 (0)4532 / 400 49-10

Email: info@mk-dent.com

www.MK-dent.com