

# DISCOVER

**MK-dent**<sup>®</sup>  
GERMANY



## PRIME LINE

HIGH & LOW SPEED HANDPIECES

[www.MK-dent.com](http://www.MK-dent.com)

# PRIME LINE



Max. 55dB:

## Quiet like no Other!

- Prime Grip
- Titanium coating
- Pushbutton system
- Premium spindle with carbide-insert
- Ceramic bearings
- Glass rod light conductor
- Pressure reduction system
- Non-return valves

Available in the following versions:

	HP21KL	HP22KL
Head size	Standard	Small
Light	✓	✓
Power (Watt)	22 W	18 W
Spray system	4-hole	3-hole
Connection	Multiflex®	Multiflex®

# PRIME LS LINE

• KaVo® compatible

• Titanium coating



	LPH11	LPH16	LPH21
Transmission	1:1	1:1	2:1
Label	blue	blue	green
Max. speed	40.000 rpm	40.000 rpm	20.000 rpm
Pushbutton	✓	✓	✓
Suitable burs	Contra angle burs	FG burs	Contra angle burs

• Prime Grip

• INTRAmatic® connection

• Stainless steel case with titanium coating



	LP11L	LP27L	LP74L
Light	✓	✓	✓
Transmission	1:1	2,7:1	7,4:1
Label	blue	green	green
Spray system	internal	internal	internal
Max. speed	40.000 rpm	14.814 rpm	5.406 rpm
Connection	INTRAmatic®	INTRAmatic®	INTRAmatic®
Suitable heads	MK-dent & KaVo® contra angle heads	MK-dent & KaVo® contra angle heads	MK-dent & KaVo® contra angle heads

*MK-dent*<sup>®</sup>  
GERMANY



*MK-dent*<sup>®</sup>  
GERMANY

Eichenweg 7b  
D-22941 Bargteheide, Germany

Phone: 0049 (0)4532 / 400 49-0

Fax: 0049 (0)4532 / 400 49-10

Email: [info@mk-dent.com](mailto:info@mk-dent.com)

[www.MK-dent.com](http://www.MK-dent.com)