

# DISCOVER

**MK-dent**<sup>®</sup>  
GERMANY



## PROPHYLAXIS

AIR FLOW HANDPIECES | PROPHY HANDPIECES | SCALER

[www.MK-dent.com](http://www.MK-dent.com)

# PROPHY LINE



## Strong Partner in Prophylaxis: Light, efficient, universal!

- Air-driven dental polishing device **PR1011**
- Available for the following brands:

**K** KaVo® **W** W&H® **S** Sirona® **NT** NSK® **B** Bien-Air®

- Ergonomically shaped, light-weight and compact
- 360° rotatable outer handle RotoGrip®-system for simple and fatigue-proof operation

### PROPHY POWDER



Prophylaxis cleaning powder for Prophy Line and other air flow handpieces, fine grain size (200 mesh)

Available flavors:	<b>PP1011</b>	Mint
	<b>PP1012</b>	Lemon
	<b>PP1013</b>	Strawberry
Content:	300g	

# SCALER

- Ultrasonic scaling device for prophylaxis
- Ergonomic shape for excellent grip
- Available for the following brands:

**E** EMS®    **SA** Satelec®

- Lean steel sleeve for a long durability
- Available with a 6-LED light ring for optimum lighting
- Painfree treatments for the patient through linear movements of the instrument



SC21E



SC21EL



SC21S



SC21SL

	SC21E	SC21EL	SC21S	SC21SL
Light	—	✓	—	✓
Frequency	28 - 32 kHz	28 - 32 kHz	28 - 32 kHz	28 - 32 kHz
Funnel cover	Steel	Steel	Steel	Steel
Connection	EMS®	EMS®	Acteon®/Satelec®	Acteon®/Satelec®

## TORQUE KEY



### SC21KEY - Torque key for insertion and storage of working tips

- Combination of torque key and key holder for neat storage and sterilization of working tips
- No contact required during handling of the instrument - therefore no risk of injury or contamination
- A defined torque of 0,3 - 0,4 Nm prevents an overwinding of the thread
- Autoclavable at up to 135 °C

# BASIC LS LINE

## CONTRA ANGLE HEADS



**LBH11**



**LBH11P**

Transmission	1:1	1:1
Max. speed	20.000 rpm	5.000 rpm
Pushbutton	✓	—
Suitable attachments	Contra angle burs	Snap-on prophy cups

## CONTRA ANGLE & PROPHY HANDPIECES

	<b>LB41/11</b>	<b>LB41/11P</b>
Transmission	4:1	4:1
Max. speed	5.000 rpm	5.000 rpm
Spray system	external clip-on	—
Label	green	green
Connection	INTRAmatic®	INTRAmatic®
Suitable attachments	Contra angle burs	Snap-on prophy cups



FOR MORE HIGH & LOW SPEED HANDPIECES VISIT OUR WEBSITE  
[www.MK-dent.com](http://www.MK-dent.com)