

DISCOVER

MK-dent[®]
GERMANY



QUICK CONNECTORS & MOTORS

www.MK-dent.com

QUICK CONNECTORS



6-hole quick connector with spray regulation system for high speed handpieces

Available:	QC6016KW	with LED
	QC5016KW	with XENON
Connection:	MK-dent & KaVo Multiflex®	



6-hole quick connector for high speed handpieces

Available:	QC6016K	with LED
	QC5016K	with XENON
Connection:	MK-dent & KaVo Multiflex®	



4-hole quick connector for high speed handpieces

Available:	QC4014K	without light
Connection:	MK-dent & KaVo Multiflex®	



6-hole quick connector with spray regulation system for high speed handpieces

Available:	QC6016S	with LED
	QC5016S	with XENON
Connection:	Sirona®	



Quick connector for high speed handpieces

Available:	QC5016W	6-hole with power supply
	QC4014W	4-hole without power supply
Connection:	W&H Roto Quick®	



Quick connector for high speed handpieces

Available:	QC6016NT	6-hole with LED
	QC5016NT	6-hole with XENON
	QC4014NT	4-hole without light
Connection:	NSK MachLite®	

For all Straight & Contra Angle Handpieces:
Powerful E-drive!



- Electric motor EM1190K
- Max. 40.000 rpm
- LED
- Spray regulation system
- Guided interior separate coolants
- Compatible with KaVo® 190, 196 and K200 connections
- Short INTRAmatic® 2 connection, suitable for all common instruments

AIRMOTORS

AIRMOTOR WITH INTERNAL WATER SUPPLY



AM1004 Airmotor with internal water supply

- Max. 20.000 rpm
- Infinitely adjustable right and left drive
- Guided interior separate coolants
- Standard INTRAmatic® connection

Connection: 4-hole

AIRMOTORS WITH EXTERNAL WATER SUPPLY



AM0012 Airmotor with external water supply

- Max. 20.000 rpm
- Standard INTRAmatic® connection

Connection: 2-hole



AM0014 Airmotor with external water supply

- Max. 20.000 rpm
- Standard INTRAmatic® connection

Connection: 4-hole

AIRMOTORS WITHOUT WATER SUPPLY



AM0002 Airmotor without water supply

- Max. 20.000 rpm
- Standard INTRAmatic® connection

Connection: 2-hole



AM0004 Airmotor without water supply

- Max. 20.000 rpm
- Standard INTRAmatic® connection

Connection: 4-hole